



Cannuflow ClearVu® Product Update Notice 12/18 for REF CV-27S and CV-27B


You will notice important user-requested upgrades to ClearVu cannulas you will be receiving, starting in 2019. The ordering reference number will not change.

New features:

- 1) Attractive new design with flat ergonomically designed grip area on the valve body for positive control when inserting the cannula
- 2) Proximal ribs to keep the ClearVu from backing out during manipulation of the knee during surgery
- 3) Longer luer connector to accommodate a wider range of tube set fittings such as luer connectors and funnel connectors
- 4) Positive valve stop for accurate fluid flow control
 - Turn right/counterclockwise to open flow valve
 - Turn left/clockwise to close

Other than these improvements the ClearVu is exactly the same product as the flexible, atraumatic inflow/outflow cannula you have been using and have come to rely on for superior fluid flow in knee and joint surgery. ClearVu is still the only flexible, atraumatic cannula progressively flexible available today that is clinically proven to not scratch and damage cartilage, or cause the post operative 3rd portal pain that metal cannulas can cause. The ClearVu is still the only cannula available that is low profile, and provides abundant fluid flow (> 1 liter per minute) and does not kink or clog.

Cannuflow provides free samples for qualified product trials.



Cannuflow Inc. 489 Division Street Campbell CA 95008
Phone 408-764-0220 fax 408-280-0323